



ATP

Colorimetric Microplate Assay Kit

User Manual

Catalog # CAK1120

(Version 1.5B)

Detection and Quantification of ATP content in Urine, Serum, Plasma,
Tissue extracts, Cell lysate, Cell culture media and Other biological
fluids Samples.

For research use only. Not for diagnostic or therapeutic procedures.

I. INTRODUCTION.....	2
II. KIT COMPONENTS.....	3
III. MATERIALS REQUIRED BUT NOT PROVIDED.....	4
IV. REAGENT PREPARATION.....	5
V. SAMPLE PREPARATION.....	6
VI. ASSAY PROCEDURE.....	7
VII. CALCULATION.....	8
VIII. TYPICAL DATA.....	9

I. INTRODUCTION

ATP (Adenosine 5'-triphosphate) is the chemical energy for cellular metabolism and is often referred to as "energy currency" of the cell. ATP is produced only in living cells during photosynthesis and cellular respiration and consumed in cellular processes including biosynthetic reactions, motility and cell division. It is a key indicator of cellular activity and has been utilized as a measure of cell viability and cytotoxicity in research and drug discovery.

ATP Colorimetric Microplate Assay Kit is a sensitive assay for determining ATP in various samples. ATP concentration is determined by creatine kinase and creatine. The reaction products can be measured at a colorimetric readout at 660 nm.

III.KIT COMPONENTS

Component	Volume	Storage
96-WellMicroplate	1 plate	
Assay Buffer	30 mlx 4	4 °C
Enzyme	Powder x 1	-20 °C, keep in dark
Reaction Buffer	1 mlx 1	4 °C
Substrate	Powder x 1	-20 °C, keep in dark
Dye Reagent I	Powder x 1	4 °C
Dye Reagent II	Powder x 1	4 °C
Dye Reagent III	10 ml x 1	4 °C
Standard	Powder x 1	-20 °C
Plate Adhesive Strips	3 Strips	
Technical Manual	1 Manual	

III. MATERIALS REQUIRED BUT NOT PROVIDED

1. Microplate reader to read absorbance at 660 nm
2. Distilled water
3. Pipettor, multi-channel pipettor
4. Pipette tips
5. Mortar
6. Centrifuge
7. Timer
8. Ice

IV. REAGENT PREPARATION

Important Recommendation: *Use disposable plastic tube to avoid phosphorus pollution.*

Standard: Briefly centrifuge prior to opening. Dissolve in 1 ml distilled water to generate 5 mmol/L of standard stock solution. Store at 4 °C for 1 week or -20 °C for 1-2 months after reconstitution. Then dilute to 2.5 mmol/L standard top solution by adding 500 µl stock solution into 500 µl distilled water. Perform 2-fold serial dilutions of the top standard solution using distilled water to make the standard curve. The concentration of standard curve could be 2.5/1.25/0.625/0.312/0.156/0.078/0.039 mmol/L.

Enzyme: Briefly centrifuge prior to opening. Add 1 ml distilled water to dissolve before use. Store at 4 °C for 1-2 days or -20 °C for 2 weeks.

Substrate: Add 6 ml distilled water and heat to dissolve before use. Store at 4 °C for 1 week or -20 °C for 2-3 months after reconstitution.

Dye Reagent: Add 5 ml Dye Reagent III into Dye Reagent I and 1 ml Dye Reagent III into Dye Reagent II respectively to dissolve. Transfer all Dye Reagent II into Dye Reagent III, mix, then transfer all Dye Reagent I into Dye Reagent III (Must follow this step). The mixed Dye Reagent may store at -20 °C for 2-3 weeks. ***Note:** It should be yellow. If colorless, the solution is failure. If blue, the solution is polluted. This solution should be prepared before use.

Note: Divide into small aliquots to avoid repeated freeze-thaw cycles.

V. SAMPLE PREPARATION

1. For cell and bacteria samples

Collect cell or bacteria into centrifuge tube, discard the supernatant after centrifugation, add 1ml Assay Buffer for 5×10^6 cell or bacteria, sonicate (with power 20%, sonicate 3s, interval 10s, repeat 30 times), centrifuged at 8000g 4°C for 10 minutes, take the supernatant into a new centrifuge tube, and keep it on ice for detection.

2. For tissue samples

Weigh out 0.1 g tissue, homogenize with 1ml Assay Buffer, sonicate (with power 20%, sonicate 3s, interval 10s, repeat 30 times), centrifuged at 8000g 4°C for 10 minutes, take the supernatant into a new centrifuge tube, and keep it on ice for detection.

3. For serum or plasma samples

Detect directly.

VI. ASSAY PROCEDURE

Warm all the reagents to 37°C before use.

Add following reagents into the microplate:

Reagent*	Sample**	Standard	Blank
Sample	20 µl	--	--
Standard	--	20 µl	--
Distilled water	--	--	20 µl
Substrate	60 µl	60 µl	60 µl
Reaction Buffer	10 µl	10 µl	10 µl
Enzyme	10 µl	10 µl	10 µl
Mix, put it in the oven, 37°C for 30 minutes.			
Dye Reagent Working Solution	100 µl	100 µl	100 µl
Mix, room temperature for 20 minutes, record absorbance measured at 660 nm.			

Note:

*Reagents must be added sequentially and should not be premixed prior to addition.

**The concentrations can vary over a wide range depending on the different samples.

For unknown samples, we recommend doing a pilot experiment & testing several doses to ensure the readings are within the standard curve range.

**For colored samples, we recommend setting a parallel sample background control well with same volume of sample only. Other reagents were replaced by distilled water to the same total volume. Subtract the OD value of the sample background control from the OD value of the sample to correct for interference from the sample's own color.

VII. CALCULATION

1. Calculate the sample concentration in ASSAY PROCEDURE according to the slope of the standard curve

$$C = \frac{(OD_{\text{Sample}} - OD_{\text{Blank}}) - \text{Intercept}}{\text{Slope}} \times n (\text{mmol/L})$$

Calculate the initial concentration according to sample preparation procedure.

2. According to one point of the standard OD and concentration

2.1 According to the protein concentration of sample

$$C = \frac{(C_{\text{Standard}} \times V_{\text{Standard}}) \times (OD_{\text{Sample}} - OD_{\text{Blank}})}{(OD_{\text{Standard}} - OD_{\text{Blank}}) \times C_{\text{Protein}} \times V_{\text{Sample}}} (\mu\text{mol/mg})$$

2.2 According to the quantity of cells or bacteria

$$C = \frac{(C_{\text{Standard}} \times V_{\text{Standard}}) \times (OD_{\text{Sample}} - OD_{\text{Blank}})}{(OD_{\text{Standard}} - OD_{\text{Blank}}) \times N \times (V_{\text{Sample}} / V_{\text{Assay}})} (\mu\text{mol}/10^4)$$

2.3 According to the weight of sample

$$C = \frac{(C_{\text{Standard}} \times V_{\text{Standard}}) \times (OD_{\text{Sample}} - OD_{\text{Blank}})}{(OD_{\text{Standard}} - OD_{\text{Blank}}) \times W \times (V_{\text{Sample}} / V_{\text{Assay}})} (\mu\text{mol/g})$$

2.4 According to the volume of sample

$$C = \frac{(C_{\text{Standard}} \times V_{\text{Standard}}) \times (OD_{\text{Sample}} - OD_{\text{Blank}})}{(OD_{\text{Standard}} - OD_{\text{Blank}}) \times V_{\text{Sample}}} (\mu\text{mol/ml})$$

Slope: the absorbance slope of standard curve

n: the dilution factor

C_{Standard} : the standard concentration, mmol/L = $\mu\text{mol/ml}$

V_{Standard} : the volume of standard in assay procedure, μl

V_{Sample} : the volume of sample in assay procedure, μl

V_{Assay} : the volume of Assay Buffer, μl


C_{Protein} : the sample protein concentration, mg/ml

W: the weight of sample, g

N: the quantity of cell or bacteria, 10^4

VIII. TYPICAL DATA

The standard curve is for demonstration only. A standard curve must be run with each assay.



Detection Range: 0.02mmol/L -2.5mmol/L