

## Product Data Sheet


## Human PEX10 Knockdown HEK293T Cell Lysate

| Catalog #                 | Source                                                                                                                                                                                                                                                                                                                                                                                                                                             | Reactivity | Applications |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------------|
| CKD2105                   | HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                            | Human      |              |
| <b>Description</b>        | Human PEX10 knockdown HEK293T cell lysate is engineered by our optimized transduction of the specific shRNA.                                                                                                                                                                                                                                                                                                                                       |            |              |
| <b>Form</b>               | Liquid in 50 mM Tris-HCl, 2% SDS, 10% glycerol, 100 mM DTT, 0.01% Bromophenol blue                                                                                                                                                                                                                                                                                                                                                                 |            |              |
| <b>Gene Symbol</b>        | PEX10                                                                                                                                                                                                                                                                                                                                                                                                                                              |            |              |
| <b>Alternative Names</b>  | RNF69; Peroxisome biogenesis factor 10; Peroxin-10; Peroxisomal biogenesis factor 10; Peroxisome assembly protein 10; RING finger protein 69                                                                                                                                                                                                                                                                                                       |            |              |
| <b>Entrez Gene</b>        | 5192                                                                                                                                                                                                                                                                                                                                                                                                                                               |            |              |
| <b>SwissProt</b>          | O60683                                                                                                                                                                                                                                                                                                                                                                                                                                             |            |              |
| <b>KO/KD Validation</b>   | RT-qPCR                                                                                                                                                                                                                                                                                                                                                                                                                                            |            |              |
| <b>Cell Line</b>          | HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                            |            |              |
| <b>Components</b>         | Human PEX10 knockdown HEK293T cell lysate x 1 vial; shControl HEK293T cell lysate x 1 vial                                                                                                                                                                                                                                                                                                                                                         |            |              |
| <b>Directions for Use</b> | Thaw the vial on ice or at room temperature. Mix gently by pipetting or brief vortexing. Briefly centrifuge to collect the sample. Load the desired amount of lysate directly onto an SDS-PAGE gel. Perform SDS-PAGE and Western blot analysis using standard laboratory protocols. Note: This lysate has been pre-denatured in SDS sample buffer containing DTT and does not require additional reducing agent or boiling before electrophoresis. |            |              |
| <b>Storage/Stability</b>  | Shipped at 4°C. Upon delivery aliquot and store at -80°C for one year. Avoid freeze/thaw cycles.                                                                                                                                                                                                                                                                                                                                                   |            |              |


## COHESION BIOSCIENCES LIMITED

WEB  
[www.cohesionbio.com](http://www.cohesionbio.com)ORDER  
[order@cohesionbio.com](mailto:order@cohesionbio.com)SUPPORT  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)CUSTOM  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)


# Product Data Sheet




RT-qPCR analysis of PEX10 knockdown efficiency. Compared with the 293T shRNA control (left), PEX10 mRNA expression in the stable shRNA knockdown cell line (right) decreased, and the knockdown efficiency was 75%.



RT-qPCR amplification curves. The Ct values of the knockout sample significantly increased, confirming the effective knockdown of PEX10 mRNA. [Red curve (knockdown sample), Blue curve (shRNA control), Purple and Green curve (reference gene GAPDH)]



RT-qPCR melt curve analysis demonstrated highly specific PCR amplification product without detectable nonspecific amplification or primer-dimer formation. [Red curve (knockdown sample), Blue curve (shRNA control), Purple and Green curve (reference gene GAPDH)]



This picture presents the raw fluorescence data (relative fluorescence units, RFU) used to generate the melt curves. [Red curve (knockdown sample), Blue curve (shRNA control), Purple and Green curve (reference gene GAPDH)]

**COHESION BIOSCIENCES LIMITED**