


## Product Data Sheet

## Human SND1 Knockdown HEK293T Cell Lysate


| Catalog #                 | Source                                                                                                                                                                                                                                                                                                                                                                                                                                             | Reactivity | Applications |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------------|
| CKD2121                   | HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                            | Human      |              |
| <b>Description</b>        | Human SND1 knockdown HEK293T cell lysate is engineered by our optimized transduction of the specific shRNA.                                                                                                                                                                                                                                                                                                                                        |            |              |
| <b>Form</b>               | Liquid in 50 mM Tris-HCl, 2% SDS, 10% glycerol, 100 mM DTT, 0.01% Bromophenol blue                                                                                                                                                                                                                                                                                                                                                                 |            |              |
| <b>Gene Symbol</b>        | SND1                                                                                                                                                                                                                                                                                                                                                                                                                                               |            |              |
| <b>Alternative Names</b>  | TDRD11; Staphylococcal nuclease domain-containing protein 1; 100 kDa coactivator; EBNA2 coactivator p100; Tudor domain-containing protein 11; p100 co-activator                                                                                                                                                                                                                                                                                    |            |              |
| <b>Entrez Gene</b>        | 27044                                                                                                                                                                                                                                                                                                                                                                                                                                              |            |              |
| <b>SwissProt</b>          | Q7KZF4                                                                                                                                                                                                                                                                                                                                                                                                                                             |            |              |
| <b>KO/KD Validation</b>   | RT-qPCR                                                                                                                                                                                                                                                                                                                                                                                                                                            |            |              |
| <b>Cell Line</b>          | HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                            |            |              |
| <b>Components</b>         | Human SND1 knockdown HEK293T cell lysate x 1 vial; shControl HEK293T cell lysate x 1 vial                                                                                                                                                                                                                                                                                                                                                          |            |              |
| <b>Directions for Use</b> | Thaw the vial on ice or at room temperature. Mix gently by pipetting or brief vortexing. Briefly centrifuge to collect the sample. Load the desired amount of lysate directly onto an SDS-PAGE gel. Perform SDS-PAGE and Western blot analysis using standard laboratory protocols. Note: This lysate has been pre-denatured in SDS sample buffer containing DTT and does not require additional reducing agent or boiling before electrophoresis. |            |              |
| <b>Storage/Stability</b>  | Shipped at 4°C. Upon delivery aliquot and store at -80°C for one year. Avoid freeze/thaw cycles.                                                                                                                                                                                                                                                                                                                                                   |            |              |

## COHESION BIOSCIENCES LIMITED


# Product Data Sheet




RT-qPCR analysis of SND1 knockdown efficiency. Compared with the 293T shRNA control (left), SND1 mRNA expression in the stable shRNA knockdown cell line (right) decreased, and the knockdown efficiency was 85%.



RT-qPCR amplification curves. The Ct values of the knockout sample significantly increased, confirming the effective knockdown of SND1 mRNA. [Red curve (knockdown sample), Blue curve (shRNA control), Purple and Green curve (reference gene GAPDH)]



RT-qPCR melt curve analysis demonstrated highly specific PCR amplification product without detectable nonspecific amplification or primer-dimer formation. [Red curve (knockdown sample), Blue curve (shRNA control), Purple and Green curve (reference gene GAPDH)]



This picture presents the raw fluorescence data (relative fluorescence units, RFU) used to generate the melt curves. [Red curve (knockdown sample), Blue curve (shRNA control), Purple and Green curve (reference gene GAPDH)]

**COHESION BIOSCIENCES LIMITED**