

Product Data Sheet

Primary Antibody Diluent

Catalog #	Source	Reactivity	Applications
CRG1027		N/A	
Description	Used to dilute primary antibodies		
Form	Liquid		
Components	1% BSA, TBS, 0.05% Proclin300, pH 7.4		
Application	Primary Antibody Diluent can be used to dilute antibodies for use in WB, IHC, IF procedures.		

Storage/Stability

1. Shipping by ice packs. Store at 2-8 °C for short-term use; store at -20 °C for long-term storage. Product is stable for one year from the date of receipt. Avoid repeated freeze-thaw cycles; aliquot into single-use portions before storage.
2. After opening the product, it should be stored at 2-8 °C and use within two weeks. For longer storage, aliquot sterilely in a laminar-flow hood immediately after opening, freeze the aliquots at -20 °C. After thawing an aliquot, keep at 2-8 °C and use within two weeks; do not re-freeze.
3. Dilute antibody immediately before use. Do not store or re-use diluted antibody solutions.

COHESION BIOSCIENCES LIMITED