

Product Data Sheet

Recombinant Human FKBP7 Protein

Catalog #	Source	Reactivity	Applications
CRP1570	Human cells	H	E, WB, SDS-PAGE, MS
Description	Recombinant Human FKBP7 Protein is produced by mammalian expression system and the target gene encoding Gln24-Leu222 is expressed with a 6His tag at the C-terminus.		
Form	Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM CaCl2, 10% Glycerol, pH 7.5.		
Gene Symbol	FKBP7		
Alternative Names	FKBP23; Peptidyl-prolyl cis-trans isomerase FKBP7; PPIase FKBP7; 23 kDa FK506-binding protein; 23 kDa FKBP; FKBP-23; FK506-binding protein 7; FKBP-7; Rotamase		
Entrez Gene	51661 (Human)		
SwissProt	Q9Y680 (Human)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	QRQKKEESTE EVKIEVLHRP ENCSKTSKKG DLLNAHYDGY LAKDGSKFYC SRTQNEGHPK WFLVGVGQVI KGLDIAMTDM CPGEKRKQVI PPSFAYGKEG YAEGKIPDA TLIFEIELYA VTKGPRSIET FKQIDMDNDR QLSKAEINLY LQREFEKDEK PRDKSYQDAV LEDIFKKNDH DGDGFISPKE YNVYQHDELV DHHHHHH		
Quality Control	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be		
Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference			
Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish			

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com