

Product Data Sheet

Recombinant Human METAP1D Protein

Catalog #	Source	Reactivity	Applications
CRP1954	E. coli	H	E, WB, SDS-PAGE, MS
Description	Recombinant Human METAP1D Protein is produced by E. coli expression system and the target gene encoding Arg44-Ala335 is expressed with a 6His tag at the N-terminus, 6His tag at the C-terminus.		
Form	Liquid in a 0.2 μM filtered solution of 20mM Tris HCl, pH8.0.		
Gene Symbol	METAP1D		
Alternative Names	MAP1D; Methionine aminopeptidase 1D mitochondrial; MAP 1D; MetAP 1D; Methionyl aminopeptidase type 1D mitochondrial; Peptidase M 1D		
Entrez Gene	254042 (Human)		
SwissProt	Q6UB28 (Human)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	MGSSHHHHHH SSGLVPRGSH MRQRDISHSI VLPAAVSSAH PVPKHIKKPD YVTTGIVPDW GDSIEVKNEQ IQQLHQACQ LARHVLLLAG KSLKVDMTTE EIDALVHREI ISHNAYPSPL GYGGFPKSVQ TSVNNVLCHG IPDSRPLQDG DIINIDVTYV YNGYHGDTSQ TFLVGNVDEC GKKLVEVARR CRDEAIAACR AGAPFSVIGN TISHITHQNG FQVCPHFVGH GIGSYFHGHP EIWHHANDSD LPMEEGMAFT IEPITEGSP EFKVLEDAWT VVSLDNQRSA QFEHTVLITS RGAQILTKLP HEALEHHHHH H		
Quality Control	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be		
Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference			
Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish			

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com