

# Product Data Sheet

## Recombinant Mouse NT5E Protein

| Catalog #                                                                                                                                                                                                                                                                                                                             | Source                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | Reactivity | Applications        |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------------|
| CRP2031                                                                                                                                                                                                                                                                                                                               | Human cells                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | M          | E, WB, SDS-PAGE, MS |
| Description                                                                                                                                                                                                                                                                                                                           | Recombinant Mouse NT5E Protein is produced by mammalian expression system and the target gene encoding Trp29-Phe550 is expressed with a 6His tag at the C-terminus.                                                                                                                                                                                                                                                                                                                                                                                                                                                         |            |                     |
| Form                                                                                                                                                                                                                                                                                                                                  | Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 120mM NaCl, 4mM CaCl2, 20% Glycerol, pH 7.5.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |            |                     |
| Gene Symbol                                                                                                                                                                                                                                                                                                                           | NT5E                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |            |                     |
| Alternative Names                                                                                                                                                                                                                                                                                                                     | NT5; NTE; 5'-nucleotidase; 5'-NT; Ecto-5'-nucleotidase; CD73                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |            |                     |
| Entrez Gene                                                                                                                                                                                                                                                                                                                           | 23959 (Mouse)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |            |                     |
| SwissProt                                                                                                                                                                                                                                                                                                                             | Q61503 (Mouse)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |            |                     |
| Purity                                                                                                                                                                                                                                                                                                                                | Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |            |                     |
| Chemical Structure                                                                                                                                                                                                                                                                                                                    | WELTILHTND VHSRLEQTSD DSTKCLNASL CVGGVARLFT KVQQIRKEEP NVLFLDAGDQ<br>YQGTIWFTVY KGLEVAHFMN ILGYDAMALG NHEFDNGVEG LIDPLLRNVK FPILSANIKA<br>RGPLAHQISG LFLPSKVLSV GGEVVGIVGY TSKETPFLSN PGTNLVFEDE ISALQPEVDK<br>LKTLNVNKII ALGHSGFEMD KLIAQKVRGV DIVVGGHSNT FLYTGNPPSK EVPAGKYPII<br>VTADDGRQVP VVAYAFGKYL GYLKVEFDDK GNVITSYGNP ILLNSSIPED ATIKADINQW<br>RIKLDNYSTQ ELGRTIVYLD GSTQTCRFRE CNMGNLICDA MINNNLRHPD EMFWNHVSMC<br>IVNGGGIRSP IDEKNNGTIT WENLAAVLPF GGTFDLVQLK GSTLKKAFEH SVHRYGQSTG<br>EFLQVGGIHV VYDINRKPWN RVVQLEVLCT KCRVPIYEPL EMDKVYKVTL PSYLANGGDG<br>FQMIKDELLK HDSGDQDISV VSEYISKMKV VYPAVEGRI KFHHHHHH |            |                     |
| Quality Control                                                                                                                                                                                                                                                                                                                       | Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |            |                     |
| Directions for Use                                                                                                                                                                                                                                                                                                                    | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |            |                     |
| Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |            |                     |
| Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish                                                                                                                                                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |            |                     |

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com

## Product Data Sheet

recommended to reconstitute to a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

### Storage/Stability

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

#### WEB

[www.cohesionbio.com](http://www.cohesionbio.com)

#### ORDER

[order@cohesionbio.com](mailto:order@cohesionbio.com)

#### SUPPORT

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

#### CUSTOM

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)