

# Product Data Sheet

## Recombinant Human PMVK Protein

| Catalog #                 | Source                                                                                                                                                                                                                                                | Reactivity | Applications        |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------------|
| CRP2124                   | E. coli                                                                                                                                                                                                                                               | H          | E, WB, SDS-PAGE, MS |
| <b>Description</b>        | Recombinant Human PMVK Protein is produced by E.coli expression system and the target gene encoding Met1-Leu192 is expressed with a 6His tag at the N-terminus.                                                                                       |            |                     |
| <b>Form</b>               | Liquid in a 0.2 µM filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 7.5.                                                                                                                                                     |            |                     |
| <b>Gene Symbol</b>        | PMVK                                                                                                                                                                                                                                                  |            |                     |
| <b>Alternative Names</b>  | PMKI; Phosphomevalonate kinase; PMKase; hPMK                                                                                                                                                                                                          |            |                     |
| <b>Entrez Gene</b>        | 10654 (Human)                                                                                                                                                                                                                                         |            |                     |
| <b>SwissProt</b>          | Q15126 (Human)                                                                                                                                                                                                                                        |            |                     |
| <b>Purity</b>             | Greater than 95% as determined by reducing SDS-PAGE.                                                                                                                                                                                                  |            |                     |
| <b>Chemical Structure</b> | MGSSHHHHHH SSGLVPRGSH MAPLGGAAPRL VLLFSGKRKS GKDFVTEALQ SRLGADVCAV<br>LRLSGPLKEQ YAQEHGLNFQ RLLDTSTYKE AFRKDMIRWG EEKRQADPGF FCRKIVEGIS<br>QPIWLVS DTR RVSDIQWFRE AYGAVTQTVR VVALEQSRQQ RGWVFTPGVD DAESECGLDN<br>FGDFDWVIEN HGVEQRLEE Q LENLIEFIRS RL |            |                     |
| <b>Quality Control</b>    | Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test.                                                                                                                                                                                  |            |                     |
| <b>Directions for Use</b> | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.                    |            |                     |
| <b>Storage/Stability</b>  | Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.                                                                                                      |            |                     |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com