

## Recombinant Human PSMD10 Protein

| Catalog #                 | Source                                                                                                                                                                                                                                                            | Reactivity | Applications        |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------------|
| CRP2175                   | E. coli                                                                                                                                                                                                                                                           | H          | E, WB, SDS-PAGE, MS |
| <b>Description</b>        | Recombinant Human PSMD10 Protein is produced by E. coli expression system and the target gene encoding Met1-Gly226 is expressed.                                                                                                                                  |            |                     |
| <b>Form</b>               | Liquid in a 0.2 µM filtered solution of PBS, 10% Glycerol, pH7.4.                                                                                                                                                                                                 |            |                     |
| <b>Gene Symbol</b>        | PSMD10                                                                                                                                                                                                                                                            |            |                     |
| <b>Alternative Names</b>  | 26S proteasome non-ATPase regulatory subunit 10; 26S proteasome regulatory subunit p28; Gankyrin; p28(GANK                                                                                                                                                        |            |                     |
| <b>Entrez Gene</b>        | 5716 (Human)                                                                                                                                                                                                                                                      |            |                     |
| <b>SwissProt</b>          | O75832 (Human)                                                                                                                                                                                                                                                    |            |                     |
| <b>Purity</b>             | Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.                                                                                                                                                                                                 |            |                     |
| <b>Chemical Structure</b> | MEGCVSNLMV CNLAYSQKLE ELKESILADK SLATRTDQDS RTALHWACSA GHTEIVEFL<br>QLGVPVNDKD DAGWSPLHIA ASAGRDEIVK ALLGKGAQVN AVNQNGCTPL<br>HYAASKNRHE IAVMLLEGGA NPDAKDHYEA TAMHRAAAKG NLKMIHILLY YKASTNIQDT<br>EGNTPLHLAC DEERVEEAKL LVSQGASIYI ENKEEKTPLQ VAKGGGLGLIL KRMVEG |            |                     |
| <b>Quality Control</b>    | Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test.                                                                                                                                                                                              |            |                     |
| <b>Directions for Use</b> | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.                                |            |                     |
| <b>Storage/Stability</b>  | Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.                                                                                                                  |            |                     |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com