

# Product Data Sheet

## Recombinant Human UBE2G2 Protein

| Catalog #                 | Source                                                                                                                                                                                                                                                                                                                                                                                                                                                             | Reactivity | Applications        |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------------|
| CRP2480                   | E. coli                                                                                                                                                                                                                                                                                                                                                                                                                                                            | H          | E, WB, SDS-PAGE, MS |
| <b>Description</b>        | Recombinant Human UBE2G2 Protein is produced by E. coli expression system and the target gene encoding Met1-Leu165 is expressed with a GST tag at the N-terminus.                                                                                                                                                                                                                                                                                                  |            |                     |
| <b>Form</b>               | Liquid in a 0.2 µM filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.                                                                                                                                                                                                                                                                                                                                                                     |            |                     |
| <b>Gene Symbol</b>        | UBE2G2                                                                                                                                                                                                                                                                                                                                                                                                                                                             |            |                     |
| <b>Alternative Names</b>  | Ubiquitin-conjugating enzyme E2 G2; Ubiquitin carrier protein G2; Ubiquitin-protein ligase G2                                                                                                                                                                                                                                                                                                                                                                      |            |                     |
| <b>Entrez Gene</b>        | 7327 (Human)                                                                                                                                                                                                                                                                                                                                                                                                                                                       |            |                     |
| <b>SwissProt</b>          | P60604 (Human)                                                                                                                                                                                                                                                                                                                                                                                                                                                     |            |                     |
| <b>Purity</b>             | Greater than 95% as determined by reducing SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                                                               |            |                     |
| <b>Chemical Structure</b> | MSPILGYWKI KGLVQPTRL LLEYLEEKYEE HLYERDEGDK WRNKKFELGL EFPNLPYYID<br>GDVKLTQSMA IIRYIADKHN MLGGCPKERA EISMLEGAVL DIRYGVSRIA YSKDFETLKV<br>DFLSKLPEML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD VVLYMDPMCL DAFPKLVCFK<br>KRIEAIPIQID KYLKSSKYIA WPLQGWQATF GGGDHPPKSD LVPRGSHMAG TALKRLMAEY<br>KQLTLNPPEG IVAGPMNEEN FFEWEALIMG PEDTCFEFGV FPAILSFPLD YPLSPPKMRF<br>TCEMFHPNIY PDGRVCISIL HAPGDDPMGY ESSAERWSPV QSVEKILLSV VSMLAEPNDE<br>SGANVDASKM WRDDREQFYK IAKQIVQKSL GL |            |                     |
| <b>Quality Control</b>    | Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test.                                                                                                                                                                                                                                                                                                                                                                                               |            |                     |
| <b>Directions for Use</b> | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Please aliquot                                                                                                                                                                                                                                                                                            |            |                     |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com

## Product Data Sheet

the reconstituted solution to minimize freeze-thaw cycles.

### Storage/Stability

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

#### WEB

[www.cohesionbio.com](http://www.cohesionbio.com)

#### ORDER

[order@cohesionbio.com](mailto:order@cohesionbio.com)

#### SUPPORT

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

#### CUSTOM

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)