

Product Data Sheet

Recombinant Mouse AGR3 Protein

Catalog #	Source	Reactivity	Applications
CRP2846	Human cells	M	E, WB, SDS-PAGE, MS
Description	Recombinant Mouse AGR3 Protein is produced by mammalian expression system and the target gene encoding Ile21-Leu165 is expressed with a 6His tag at the C-terminus.		
Form	Liquid in a 0.2 µm filtered solution of 20mM Tris, 10%Glycerol, 100mM NaCl, pH8.0.		
Gene Symbol	AGR3		
Alternative Names	BCMP11; Anterior gradient protein 3 homolog; AG-3; AG3; hAG-3; Breast cancer membrane protein 11		
Entrez Gene	403205 (Mouse)		
SwissProt	Q8R3W7 (Mouse)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	IAIKKEKRPP QTLRSGWGDD ITWVQTYEEG LFHARKSNKP LMVIHHLEDC QYCQALKKEF AKNEEIQEMA QNDFIMLNLM HETTDKNLSP DGQYVPRIMF VDPSLTVRAD ITGRYSNRLY TYEPQDLPLML VDNMKKALRL IQSELVDHHH HHH		
Quality Control	Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com