

Product Data Sheet

Recombinant Mouse IDO1 Protein

Catalog #	Source	Reactivity	Applications
CRP2977	E. coli	M	E, WB, SDS-PAGE, MS
Description	Recombinant Mouse IDO1 Protein is produced by E.coli expression system and the target gene encoding Met1-Pro407 is expressed with a 6His tag at the N-terminus.		
Form	Liquid in a 0.2 µm filtered solution of 20mM Sodium Acetate, 150mM NaCl and 20% Glycerol, pH4.5.		
Gene Symbol	IDO1		
Alternative Names	IDO; INDO; Indoleamine 2,3-dioxygenase 1; IDO-1; Indoleamine-pyrrole 2,3-dioxygenase		
Entrez Gene	15930 (Mouse)		
SwissProt	P28776 (Mouse)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	MNHKVVHHHH HMAISKISPT EGSRRILEDH HIDEVDVGFAL PHPLVELPDA YSPWVLVARN LPVLIENGQL REEVEKLPTL STDGLRGHRL QRLAHLALGY ITMAYVWNRG DDDVRKVLPR NIAVPYCELS EKLGLPPILS YADCVLANWK KKDPNGPMTY ENMDILFSFP GGDCDKGFFL VSLLEIAAS PAIKAIPTVS SAVERQDLKA LEKALHDIAT SLEKAKEIFK RMRDFVDPDT FFHVLRIYLS GWKCSSLPE GLLYEGVWDT PKMFSGGSAG QSSIFQSLDV LLGIKHEAGK ESPAEFLQEM REYMPPAHRN FLFFLESAPP VREFVISRHN EDLTAKAYNEC VNGLVSVRK HLAIVDTYIM KPSKKKPTDG DKSEEPSNVE SRGTGGTNPM TFLRSVKDTT EKALLSWP		
Quality Control	Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

Storage/Stability

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com