

# Product Data Sheet

## Anti-GST-tag Antibody - Cyanine5.5

| Catalog #                | Source                                                                                                          | Reactivity | Applications |
|--------------------------|-----------------------------------------------------------------------------------------------------------------|------------|--------------|
| CSA5521                  | Mouse                                                                                                           | N/A        | IF, FC       |
| <b>Description</b>       | Mouse monoclonal antibody to GST-tag Cyanine5.5 labeled                                                         |            |              |
| <b>Immunogen</b>         | Recombinant protein corresponding to GST-tag.                                                                   |            |              |
| <b>Purification</b>      | The antibody was purified by Protein A/G affinity chromatography.                                               |            |              |
| <b>Specificity</b>       | Recognizes GST and GST tag fusion proteins.                                                                     |            |              |
| <b>Clonality</b>         | Monoclonal                                                                                                      |            |              |
| <b>Conjugation</b>       | Cyanine5.5                                                                                                      |            |              |
| <b>Form</b>              | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide |            |              |
| <b>Dilution</b>          | IF (1/500 - 1/2000), FC (1/1000 - 1/4000)                                                                       |            |              |
| <b>Fluorophore</b>       | Ex/Em = 674/690 nm                                                                                              |            |              |
| <b>Storage/Stability</b> | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                |            |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)