

Product Data Sheet

Anti-FLAG-tag Antibody - AcalephFluor680

Catalog #	Source	Reactivity	Applications
CSA5711	Rabbit	N/A	IF, FC
Description	Rabbit polyclonal antibody to FLAG-tag AcalephFluor680 labeled		
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of FLAG-tag.		
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.		
Specificity	Recognizes C-terminal, internal, and N-terminal FLAG-tag fusion proteins.		
Clonality	Polyclonal		
Conjugation	AcalephFluor680		
Form	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide		
Dilution	IF (1/500 - 1/2000), FC (1/1000 - 1/4000)		
Fluorophore	Ex/Em = 680/701 nm		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com