

Product Data Sheet

Homocysteine Lyase

Catalog #	Source	Reactivity	Applications
CZM1009	Microorganism		
Description	Recombinant Homocysteine Lyase from Microorganism		
Specificity	>30 U/mg		
Form	Lyophilized powder		
Alternative Names	EC 4.4.1.2; Hcyase; L-homocysteine hydrogen-sulfide-lyase (deaminating 2-oxobutanoate-forming); homocysteine desulfurase; L-homocysteine hydrogen-sulfide-lyase; Homocysteine desulfhydrase; Homocysteine Lyase		
CAS Number	9024-41-3		
Molecular Weight			
Purity	>90%		
Unit Definition			
Storage/Stability	Shipped at 4° C. Store at -20° C for 2 years.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com