

Product Data Sheet

S-Adenosyl-L-Homocysteine Hydrolase

Catalog #	Source	Reactivity	Applications
CZM1040	Microorganism		
Description	Recombinant S-Adenosyl-L-Homocysteine Hydrolase from Microorganism		
Specificity	>1 U/mg		
Form	White Lyophilized powder or Liquid		
Alternative Names	EC 3.3.1.1; SAHH; Adenosylhomocysteinase; S-adenosylhomocysteine synthase; S-adenosylhomocysteine hydrolase; adenosylhomocysteine hydrolase (ambiguous); S-adenosylhomocysteinase; SAHase; AdoHcyase		
CAS Number	9025-54-1		
Molecular Weight			
Purity	>90%		
Unit Definition			
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com