

Product Data Sheet

Isocitrate Dehydrogenase

Catalog #	Source	Reactivity	Applications
CZM1053	E.coli		
Description	Recombinant Isocitrate Dehydrogenase from Microorganism		
Specificity	>80 U/mg protein, >40 U/mg powder		
Form	White lyophilized powder		
Alternative Names	EC 1.1.1.42; ICDH; isocitrate dehydrogenase (NADP+); oxalosuccinate decarboxylase; oxalosuccinic decarboxylase; isocitrate (NADP) dehydrogenase; isocitrate (nicotinamide adenine dinucleotide phosphate) dehydrogenase; NADP-specific isocitrate dehydrogenase;		
CAS Number	9028-48-2		
Molecular Weight			
Purity	>90%		
Unit Definition	One unit is defined as the amount of enzyme will convert 1 micromole isocitrate to alpha-ketoglutarate per min at pH 7.4, 37° C.		
Storage/Stability	Shipped at 4° C. Store at -20° C for 2 years.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com